

## CHEM SHARP™ TUNGSTEN SHARPENER



### **A safe alternative to tungsten grinding.**

Chemically mills a perfect point on pure, thoriated, ceriated, and lanthanated tungsten. Average time to re-point tungsten is approximately one minute at the work station. Compared to traditional sharpening methods, chemical re-pointing extends tungsten life by 400%.

**Description** 6 x 4.5 oz. Jars

**Unit** Case

**Features** Biodegradable, Non-Flammable, Non-Toxic

[600B SDS](#)

**SKU:** N/A | **Categories:** [Chemical Tungsten Sharpening Products](#), [Chemicals](#), [GREEN LINE](#) |



